# plasma care®

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[www.terraplasma-medical.com](http://www.terraplasma-medical.com)

CE0197

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## 12. Emission & Immunity Tests (EN)

| Phenomenon                                 | Limit / Compliance Level                | Electromagnetic environment - guidance                                                                                                                                                                                                                                      |
|--------------------------------------------|-----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted and radiated emission            | CISPR 11, Group 1, Class B              | Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                                                                                 |
| Harmonic current emissions (IEC 61000-3-2) | Class A                                 | Device is suitable for direct connection to the public low-voltage power supply network.                                                                                                                                                                                    |
| Voltage fluctuations and flickering        | —                                       |                                                                                                                                                                                                                                                                             |
| Radiated RF EM field (IEC 61000-4-3)       | 80-2700 MHz;<br>1kHz AM 80 %;<br>10 V/m | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance |

|                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                |
|---------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Radiated RF<br>EM filed<br>(IEC 61000-4-3)                                            | 80-2700 MHz;<br>1kHz AM 80 %;<br>10 V/m                                                                                                                                                                                                                                                                                                                                                                               | $d = 1.2\sqrt{P}$ for 80 MHz<br>to 800 MHz<br>$d = 2.3\sqrt{P}$ for 800 MHz<br>to 2,7 GHz<br>where P is the maximum output<br>power rating of the transmitter<br>in watts (W) according to the<br>transmitter manufacturer and d<br>is the recommended separa-<br>tion distance in meters (m). |
| Proximity fields form<br>RF wireless commu-<br>nications equipment<br>(IEC 61000-4-3) | 385 MHz; Pulse Modu-<br>lation: 18 Hz; 27 V/m<br>450 MHz, FM + 5<br>Hz deviation: 1 kHz<br>sine; 28 V/m<br>710, 745, 780 MHz;<br>Pulse Modulation:<br>217 Hz; 9 V/m<br>810, 870, 930 MHz;<br>Pulse Modulation:<br>18 Hz; 28 V/m<br>1720, 1845, 1970 MHz;<br>Pulse Modulation:<br>217 Hz; 28 V/m<br>2450 MHz; Pulse<br>Modulation: 217<br>Hz; 28 V/m;<br>5240, 5500, 5785<br>MHz; Pulse Modula-<br>tion: 217 Hz; 9 V/m | Portable and mobile RF<br>communications equipment<br>should be used no closer to<br>any part of the device, including<br>cables, than the recommended<br>separation distance<br>30 cm.                                                                                                        |

| Phenomenon                                                                 | Limit / Compliance Level                                                                                                                                                                                                                                                                                                                                                   | Electromagnetic environment - guidance                                                                                                                                |
|----------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Proximity fields from RF wireless communications equipment (IEC 61000-4-3) | 385 MHz; Pulse Modulation: 18 Hz; 27 V/m<br>450 MHz, FM + 5 Hz deviation: 1 kHz sine; 28 V/m<br>710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m<br>810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m<br>1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m<br>2450 MHz; Pulse Modulation: 217 Hz; 28 V/m;<br>5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance 30 cm. |
| Electrical fast transients / bursts (IEC 61000-4-4)                        | Power lines: 2 kV; 100 kHz repetition frequency                                                                                                                                                                                                                                                                                                                            | Mains power quality should be that of a typical environment.                                                                                                          |
| Surges (IEC 61000-4-5)                                                     | L-N: 1kV at 0°,90°,180°,270°                                                                                                                                                                                                                                                                                                                                               | Mains power quality should be that of a typical environment.                                                                                                          |

|                                                                |                                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------------------------------------------|----------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted disturbances induced by RF fields<br>(IEC 61000-4-6) | 0.15-80 MHz; 1kHz<br>AM 80 %; 3 Vrms,<br>6 Vrms in ISM and<br>amateur radio band | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br>Recommended separation distance<br>$d = 1.2\sqrt{P}$ for 150 kHz to 80MHz<br>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). |
| Rated power frequency magnetic fields<br>(IEC 61000-4-8)       | 30 A/m, 50 Hz                                                                    | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                                                                                                                                                                                                                                    |

| Phenomenon                                               | Limit / Compliance Level                                                                                                                                   | Electromagnetic environment - guidance                                                                                                                                                                                                                 |
|----------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Voltage dips / Voltage interruptions<br>(IEC 61000-4-11) | 0 % UT for 0.5 cycle at 0°,45°,90°,135°,180°, 225°,270°,315°<br>0 % UT for 1 cycle at 0°<br>70 % UT for 25/30 cycles at 0°<br>0 % UT for 250/300 cycles 0° | Mains power quality should be that of a typical environment.<br>If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or battery. |
| Proximity magnetic fields<br>(IEC 61000-4-39)            | 30 kHz; 8 A/m; CW<br>134 kHz; 65 A/m; 50 % PM with 2,1 kHz<br>13,56 MHz; 7,5 A/m;<br>50 % PM with 50 kHz                                                   | RFID, WPT and similar equipment should be used no closer to any part of the device, including cables, than the recommended separation distance 15 cm.                                                                                                  |

| Transmission- and receiving frequencies of the plasma care |                            |
|------------------------------------------------------------|----------------------------|
| Receiving frequency                                        | 13,56 MHz                  |
| Preferred frequency                                        | 13,56 MHz                  |
| Receiver bandwidth                                         | 14 kHz                     |
| Transmission frequency                                     | 13,56 MHz                  |
| Transmitter bandwidth                                      | 14 kHz                     |
| Modulation of frequency                                    | Amplitude modulation (ASK) |
| Effective radiated power                                   | 200 mW                     |

| Characteristics of the plasma care wireless charging module |               |
|-------------------------------------------------------------|---------------|
| Input Voltage                                               | 4,75 – 13 V   |
| Output Voltage                                              | 15W           |
| Frequenzy                                                   | 110 – 148 KHz |
| Min Load                                                    | 0,1A          |
| Max Load                                                    | 1,25A         |
| Peak Load                                                   | 1,5A          |
| Storage High Temp                                           | 16h at 60 °C  |
| Storage Low Temp                                            | 16h at -20 °C |
| Operation High Temp                                         | 8h at 40 °C   |
| Operation Low Temp                                          | 8h at -20 °C  |

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#### **terrapiasma medical GmbH**